



KEYLAJET®

Low-Foaming Chelating High Alkaline Liquid

Sample Trace Analysis Undiluted Concentrate

ELEMENT	DETECTED	DETECTION LIMIT
Chlorine	NONE	Less than 100 PPM
Fluorine	NONE	Less than 10 PPM
Bromine	NONE	Less than 10 PPM
Iodine	NONE	Less than 10 PPM
Boron	NONE	Less than 10 PPM
Lead	NONE	Less than 10 PPM
Zinc	NONE	Less than 10 PPM
Cadmium	NONE	Less than 50 PPM
Arsenic	NONE	Less than 50 PPM
Antimony	NONE	Less than 50 PPM
Bismuth	NONE	Less than 10 PPM
Tin	NONE	Less than 10 PPM
Sulfur	NONE	Less than 100 PPM
Mercury	NONE	Less than 0.5 PPM
Phosphorous	An ingredient in KEYLAJET present at 1.6% by weight	
Potassium	An ingredient in KEYLAJET present at 11–16% by weight	
Sodium	An ingredient in KEYLAJET present at 1–5% by weight	

Note: KEYLAJET detergent is recommended for use at 1–10% solutions.