



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 11/16/2021

CERTIFICATE OF ANALYSIS

ALCOJET

Manufactured by: Alconox Inc.

Sub lot No: J1L1 J1L1A, J1L1B, J1L1C, J1L1D

Sub lot numbers of Master lot No: J1H0

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 11/11/2021

Expires end of: 11/23

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:

OK

White powder

Solubility: (1% solution at 25 deg C,
with agitation)

COMP

Complete

pH:(1% solution in DI water)

11.6

11.0 - 12.5

Surface Tension:(1% solution Dyne/cm)

32.7

31.0 - 36.0

Foam (1% in DI water)

NIL

NIL

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.