



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 2/4/2022

CERTIFICATE OF ANALYSIS TERGAZYME

Manufactured by: Alconox Inc.

Sub lot No: G1B2 G1B2A, G1B2B, G1B2C, G1B2D

Sub lot numbers of Master lot No: G1E0

Manuf. Cat. no.: 1301, 1303, 1304, 1304-1, 1325, 1350

Date Of Manufacture: 2/2/2022

Expires end of: 02/24

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---|------|---|
| Appearance: | OK | White powder with cream colored flakes and small brown granules |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete |
| pH:(1% solution in DI water) | 9.3 | 9.0 - 10.0 |
| Surface Tension:(1% solution Dyne/cm) | 31.4 | 29.0 - 33.0 |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.