



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Factory: 61 Cornelison Ave. Jersey City, NJ 07304 USA

Date COA Created: 10/25/2022

CERTIFICATE OF ANALYSIS
CITRANOX

Manufactured by: Alconox Inc.

Sub lot No: C1k2 C1k2F, C1k2J, C1k2K, C1k2L, C1k2M

Sub lot numbers of Master lot No: C1A2

Manuf. Cat. no.: 1801, 1801-1, 1805, 1815, 1832, 1832-1, 1855

Date Of Manufacture: 10/13/2022

Expires end of: 10/24

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <u>ANALYSIS</u> | <u>RESULT</u> | <u>ACCEPTED RANGE</u> |
|--|---------------|-----------------------------|
| Appearance: | OK | Clear, olive colored liquid |
| pH: (1% solution in DI water) | 2.8 | 2.0 - 3.0 |
| Surface Tension:(1% solution Dyne/cm) | 29.8 | 29.0 - 33.0 |
| Specific Gravity: | 1.10 | 1.09 - 1.14 |

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4
residual solvent ingredients as defined by ICH Q3C*