



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 3/31/2023

CERTIFICATE OF ANALYSIS **ALCOJET**

Manufactured by: Alconox Inc.

Sub lot No: 03-29-23-11

Sub lot numbers of Master lot No: J1A2

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 3/29/2023

Expires end of: 03/25

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:

OK

White powder

Solubility: (1% solution at 25 deg C,
with agitation)

COMP

Complete

pH:(1% solution in DI water)

11.1

11.0 - 12.5

Surface Tension:(1% solution Dyne/cm)

33.7

31.0 - 36.0

Foam (1% in DI water)

NIL

NIL

Laura Falla, QA

Alconox, Inc.

***Contains no class 1, class 2, class 3, or class 4 residual
solvent ingredients as defined by ICH Q3C guidelines.***

COA NUMBER: 1042R