



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 10/6/2023

**CERTIFICATE OF ANALYSIS**  
**TERGAJET**

**Manufactured by: Alconox Inc.**

Sub lot No: 10-03-23-12

Sub lot numbers of Master lot No: T1A2

Manuf. Cat. no.: 2201, 2203, 2204, 2204-1, 2225, 2250

Date Of Manufacture: 10/3/2023

Expires end of: 10/25

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <b><u>ANALYSIS</u></b>                                | <b><u>RESULT</u></b> | <b><u>ACCEPTED RANGE</u></b> |
|---|----------------------|------------------------------|
| Appearance:   | OK                   | White powder                 |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP                 | Complete                     |
| pH:(1% solution in DI water)                          | 11.4                 | 11.0 - 12.5                  |
| Surface Tension:(1% solution Dyne/cm)                 | 31.0                 | 31.0 – 36.0                  |
| Foam (1% in DI water)                                 | NIL                  | NIL                          |

Laura Falla, QA

Alconox, Inc.

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*