



Alconox, Inc.

The Leader in Critical-Cleaning
30 Glenn St., Suite 309, White Plains, NY 10603 USA
Tel.914.948.4040 • Fax.914.948.4088
www.alconox.com • cleaning@alconox.com

Date COA Created: 11/02/2023

CERTIFICATE OF ANALYSIS
ANF019

Manufactured by: Alconox, Inc.

Sub lot No: 10-24-23-01

Sub lot numbers of Master lot No: N/A

Manuf. Cat. no.: LQNF5, LQNF55

Date Of Manufacture: 10/24/2023

Expires end of: 10/2025

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear pale yellow to amber color liquid
pH: (as is)	8.2	8.0 - 9.5
Surface Tension:(1% solution Dyne/cm)	31.8	29.0 - 33.0
Specific Gravity:	1.08	1.06 - 1.09

Laura Falla, QA

Alconox, Inc.

***Contains no class 1, class 2, class 3, or class 4 residual
solvent ingredients as defined by ICH Q3C guidelines.***